



Since 1984
Arti Surgicals

ETO Sterilizer

- 1) Fully Automatic Ethylene Oxide Sterilizer With Hmi (touch Screen Display) To Select Cycle
- 2) 100 Eto Cartridge Or Can Be Ethylene Mixed Gas
- 3) Different Type Of Pre-loaded Cycle Are Available
- 4) We Have Industrial Ethylene Oxide Sterilizer also 300ltr To 10,000 ltr
- 5) Made up of SS304
- 6) Working Temperature 121°C

ETO Sterilizer Available Sizes

Dimension	Capacity	Model
1 x 1 x 2	55 LTR	EOR112
1 x 1 x 3	85 LTR	EOR113
1 x 1 x 4	110 LTR	EOR114
1 x 1 x 4.5	125 LTR	EOR114.5
2 x 2 x 3	285 LTR	EOR223



To learn more, call us on **+91 9969572795 | +91 9768308758**
email artisurgicals@gmail.com or visit www.artisurgicals.com