



Since 1984
Arti Surgicals

Horizontal Cylindrical Autoclave



- 1) Sterilization at standard 121°C and 134°C cycles.
- 2) Set values process values of temperature pressure will be shown on display screen.
- 3) Vacuum pulsing and drying to achieve low Sterility Assurance Level (SAL).
- 4) Automatic control of temperature, pressure values and all critical parameters.
- 5) Reports auto generated in pen drive, hard disk or printer.



Horizontal Cylindrical Autoclave Available Sizes

Dimension	Working Capacity	Single Door	Double Door
16" x 24"	80 Liters	HCS1624	-
16" x 30"	100 Liters	HCS1630	HCD1630
16" x 36"	120 Liters	HCS1636	HCD1636
20" x 30"	150 Liters	HCS2030	HCD2030
20" x 36"	185 Liters	HCS2030	HCD2030
20" x 48"	248 Liters	HCS2048	HCD2048

To learn more, call us on **+91 9969572795 | +91 9768308758**
email artisurgicals@gmail.com or visit www.artisurgicals.com