



Since 1984
Arti Surgicals

Horizontal Rectangular Autoclave

- 1) Sterilization at standard 121°C and 134°C cycles.
- 2) Set values process values of temperature pressure will be shown on display screen.
- 3) Vacuum pulsing and drying to achieve low Sterility Assurance Level (SAL).
- 4) Automatic control of temperature, pressure values and all critical parameters.
- 5) Reports auto generated in pen drive, hard disk or printer.



Horizontal Rectangular Autoclave Available Sizes

Dimension	Working Capacity	Single Door	Double Door
2 x 2 x 2.25	285 Liters	HRS222.5	HRD222.5
2 x 2 x 4	450 Liters	HRS224	HRD224
2 x 2 x 5	550 Liters	HRS225	HRD225
3 x 3 x 4	1020 Liters	HRS334	HRD334



To learn more, call us on **+91 9969572795 | +91 9768308758**
email artisurgicals@gmail.com or visit www.artisurgicals.com