



Since 1984
Arti Surgicals

Vertical Autoclave Laboratory Model



Technical Data: Double Wall Construction, made up of SS304 Working Temperature 121°C @pressure 15 PSI, Foot lifting lid.

Fully Automatic Operated

- 1) High temperature Joint less silicon gasket
- 2) I.S.I. Mark SS Heater
- 3) Digital Screen Display
- 4) Settable log/sterilization time interval.
- 5) Adjustable Soak Time in 1 min to 99 min.
- 6) Low water level switch for heater protection
- 7) Automatic Steam controlling in chamber, temp maintaining, sterilization time holding, Auto-Exhaust, Cycle end.



Dimension	Working Capacity	Manual Operated	Semi-Automatic Operated	Fully Automatic Operated
10" x 18"	23 Liters	VMD1018	VSD1018	VFD1018
12" x 20"	37 Liters	VMD1220	VSD1220	VFD1220
14" x 22"	55 Liters	VMD1422	VSD1422	VFD1422
16" x 24"	80 Liters	VMD1624	VSD1624	VFD1624
18" x 24"	100 Liters	VMD1824	VSD1824	VFD1824
18" x 30"	125 Liters	VMD1830	VSD1830	VFD1830
22" x 30"	186 Liters	VMD2230	VSD2230	VFD2230
24" x 30"	222 Liters	VMD2430	VSD2430	VFD2430

To learn more, call us on **+91 9969572795 | +91 9768308758**
email artisurgicals@gmail.com or visit www.artisurgicals.com