



Since 1984
Arti Surgicals

Vertical **Autoclave** Portable Model

- 1) Single Wall Construction, made up of SS304
- 2) Working Temperature 121°C @ pressure 15 PSI.
- 3) The Gasket is made of Superior quality Silicon rubber which is easily Replaceable. Electrically operated supplied with pressure gauge, safety valve and steam release valve.



Dimension	Working Capacity	Manual Operated	Semi-Automatic Operated
12" x 12"	17 Liters	VMS1212	VSS1212
12" x 20"	37 Liters	VMS1220	VSS1220



To learn more, call us on **+91 9969572795 | +91 9768308758**
email artisurgicals@gmail.com or visit www.artisurgicals.com