



Since 1984
Arti Surgicals

Vertical Single Lever Autoclave



Top loading autoclave with single lever opening Wall conduction made up of SS304 / SS316
working Temp 121° C @ pressure 15 quick open & close lid.
Easy to operate & Time consuming locking System.



Vertical Single Lever Autoclave Available Sizes

Dimension	Working Capacity	Double Wall Cylindrical Body	Double Wall Rectangular Body	Triple Wall Cylindrical Body
12" x 20"	37 Liters	SLC1220	SLR1220	SLC1220
14" x 22"	55 Liters	SLC1422	SLR1422	-
16" x 24"	80 Liters	SLC1624	SLR1624	SLT1624
18" x 24"	100 Liters	SLC1824	SLR1824	-
18" x 30"	125 Liters	SLC1830	SLR1830	-
22" x 30"	186 Liters	SLC2230	SLR2230	-

To learn more, call us on **+91 9969572795 | +91 9768308758**
email artisurgicals@gmail.com or visit www.artisurgicals.com